# “Abnova” Capillary Blood Collection Tube (Non-Sterile)

## Instructions For Use

REF U0399

### INTENDED USE

“Abnova” Capillary Blood Collection Tube (Non-Sterile) is a capillary tube for collection and dispensing of fingertip blood for Point-of-Care (POC) tests.

“Abnova” Capillary Blood Collection Tube (Non-Sterile) are intended for *in vitro* diagnostic examination. This product should be used with Human IgM/IgG rapid Test kit for detection of human IgM and IgG antibodies in human whole blood. The product is to be used by properly trained professionals only in accordance with the instructions.

### PRINCIPLE OF USE

“Abnova” Capillary Blood Collection Tube (Non-Sterile) uses the capillary action to draw 20 µL of fingertip blood. Use a lancet (not included in this product) for fingertip skin-puncture to obtain small volume of blood. Touch the capillary tip to the blood drop to collect the blood in a horizontal position, and then press down the piston to dispense the blood on the sample zone of test strips.

### MATERIALS AND SUPPLIES PROVIDED


“Abnova” Capillary Blood Collection Tube (Non-Sterile) contains Capillary Blood Collection Tubes and Instructions for Use.

Capillary Blood Collection Tubes (200)

- 20µL capillary blood collection tube

Assay Documentation (1)

- Instructions For Use

### MATERIALS REQUIRED BUT NOT SUPPLIED

- Lancet
- Alcohol swab
- Testing strip
- Personal protective equipments

### EXPIRATION DATE

Date of manufacture, expiration date and lot number are indicated on the packaging. Expiration date is 2 years after the date of manufacture.

### STORAGE ENVIRONMENT

- Store the product at room temperature 15°C - 30°C (59°F - 86°F) and avoid direct exposure to sunlight.
- Do not freeze or store the product outside the temperature range described above.
- Do not open the product until ready to use.
- Do not use the product after the expiration date.

### WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- These instructions must be strictly followed by a trained healthcare professional to achieve blood collection and dispensing. All users have to read the instruction prior to blood collection.
- Each Capillary blood collection tube is for one-time use only. Dispose after use. DO NOT REUSE.
- DO NOT press down the piston until dispensing blood. DO NOT use the Capillary blood collection tube if the piston is loose or is pressed down before use.
- Do not touch the capillary tip directly to wounds for blood collection.
- Ensuring there are no air bubbles present during blood collection. If any bubbles appear in the capillary tube, discard it and re-collect blood with a new Capillary blood collection tube.
- Ensuring total 20µL of blood is collected. If the blood volume in the capillary tube is less than 20µL, discard it and re-collect blood with a new Capillary Blood Collection Tube.
- Freshly collected blood should be dispensing within 3 minutes. Avoid keeping the blood in the capillary tube over 3 minutes or the blood will form clots and be difficult to dispense.
- Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
- Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.
- Dispose of the potentially contaminated materials in accordance with local, regional and national regulation of biohazard.


### SPECIMEN COLLECTION

Healthcare professional need follows collection of blood specimen guidelines established by health authority or medical institutions for specimen collection.


### PROCEDURE

#### Fingertip Blood Collection


1. Clean the fingertip with an alcohol swab and prick the fingertip firmly with a lancet.

2. Hold the collection tube in a horizontal position to collect the fresh blood (DO NOT hold the collector vertically).

3. Slightly press down the piston to smoothly dispense the blood sample.

### SYMBOLS

	Do not re-use		Do not use if package is damaged		Temperature limit: 15°C - 30°C
	Consult instruction for use		Manufacturer		Catalog number
	Lot number		Use by date		Date of manufacture
	CE Mark		<i>In vitro</i> diagnostic medical device		200 tests
	Keep dry		Keep away from sunlight		Authorized Representative in the European Community

CE IVD

mdi Europa GmbH  
Langenhagener St. 71, D-30855 Langenhagen, Germany

Abnova (Taiwan) Corporation  
Zhongli Qingpu Plant No. 326-8, Sec. 4, Zhongzheng Rd.,  
Zhongli Dist., Taoyuan City 320, Taiwan (R.O.C.)

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