

Datasheet

17 alpha-hydroxyprogesterone monoclonal antibody, clone HPRO-2

Catalog Number: MAB2859

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against 17 alpha-hydroxyprogesterone.

Clone Name: HPRO-2

Immunogen: 17 alpha-hydroxyprogesterone conjugated with BSA.

Host: Mouse

Applications: ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody reacts with 17 alpha-hydroxyprogesterone-BSA conjugate and free 17 alpha;-hydroxyprogesterone. There is no cross-reactivity with BSA.

Form: Liquid

Isotype: IgG2b

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C.